




We are a team of experts driven by our urge to achieve
perfection and make a difference in the medical aesthetics
industry.



Based in Canada and branched throughout the world: mainly Spain, Cyprus and United Arab Emirates, represented as well by other agents around the globe, Celiraz International has been striving since 2015 to create and distribute the best beauty and cosmetic products using the latest technologies.



Estilad fillers are sterile medical devices, injectable, apyrogenic, absorbable, based on cross-linked hyaluronic acid derived from bacterial fermentation. The crosslinking of Hyaluronic Acid inside the product makes the filler effect more stable. Our fillers, thanks to the different concentrations of Hyaluronic Acid put inside in each product, guarantee achievement of more goals depending on each patient needs.



Estilad FINE

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments facial volume without surgery.

Estilad® Fine offers the potential to restore youthful contour, shaping and tissue fullness with a minimally invasive aesthetic procedure. Estilad® Fine is a soft tissue restorative agent and filler, it is primarily used in injection into the deep dermis and superficial layer of the subcutaneous tissue on the face to correct and remove fine to medium facial wrinkles. Estilad® Fine special formulation is suitable when treating the smile lines and helps in treating other skin problems.





Product Features

- 1 graduated sterile pre-filled syringes of 1 ml
- 1 x 1 ml and 2 needles 30 G



5 months

Treatment

- barcode
- scars from acne
- periocular wrinkles



The times indicated should be associated with the patient's lifestyle

Estilad TOUCH

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments facial volume without surgery.

Estilad® Touch is a soft tissue agent and filler, it is primarily used in injection into the deep dermis and superficial layer of the subcutaneous tissue of the face, to remove fine lines and wrinkles, correction of both cleft lip, can reduce the visible appearance of scarring also improving congenital eyelid malposition, including eyelid retraction, ectropion, euryblepharon, epiblepharon, and abnormalities associated with a shallow orbit, by injecting of Estilad® Touch in the pretarsal and/or septal regions of the affected eyelid.





Product Features

- 1 graduated sterile pre-filled syringes of 1 ml
- 1 x 1 ml and 1 needle 27 G



6 - 8 months

Treatment

- nose
- tear trough
- lips
- temples
- medium deep grooves
- medium superficial



The times indicated should be associated with the patient's lifestyle

Estilad VOLUME

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments facial volume without surgery.

Estilad® Volume is a soft tissue agent and filler, it is primarily used in injection into the deep dermis and superficial layer of the subcutaneous tissue on the face to remove deep facial wrinkles, adjust facial contours and enhance the volume of lips. Estilad® Volume offers the potential to restore youthful contour, shaping and tissue fullness with a minimally invasive aesthetic procedure.

Estilad® Volume, highly purified sodium hyaluronate crosslinked with a long chain of molecules form a coiled 3-D hydrophilic matrix network configuration that allows nutrients, oxygen and hormones to pass through freely and to protect the regular biological functions of the tissue.





Product Features

- 1 graduated sterile pre-filled syringes of 1 ml
- 1 x 1 ml and 1 needle 27 G



8 – 12 months

Treatment

- naso-genal furrows
- marionette wrinkles
- mandibular edge
- deep wrinkles
- lips large volumes
- deep furrows



The times indicated should be associated with the patient's lifestyle

Estilad VOLUMAX

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments facial volume without surgery.

Estilad® Volumax is a soft tissue restorative agent and filler, it is primarily used for injections into the deep dermis and superficial layer of the subcutaneous tissue on the face to remove facial wrinkles and adjust facial contours. Estilad® Volumax offers the potential to restore youthful contour, shaping and tissue fullness with a minimal invasive aesthetic procedure. Estilad® Volumax can correct nose deformity. By applying a small amount of filler gradually, physicians have total control over the contours of the nose and reshape it. With Estilad® Volumax, highly purified sodium hyaluronate crosslinked with a long chain of molecules form a coiled 3-D hydrophilic matrix network configuration that allows nutrients, oxygen and hormones to pass through freely and to protect the regular biological functions of the tissue.





Product Features

- 1 graduated sterile pre-filled syringes of 1 ml
- 1 x 1 ml and 1 needle 26 G



8 -12 months

Treatment

- cheekbones
- malar area
- chin
- deep furrows
- facial lipoatrophy



The times indicated should be associated with the patient's lifestyle

Estilad CONTOUR

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments body volume without surgery.

Estilad® Contour can enhance the volume and appearance of different body areas. Popular areas include the hands, décolletage, earlobes and feet. Estilad® Contour is designed for deep tissue implantation and differs in terms of the physical structure of the gel. The depth of the injection may vary from subcutaneous to supraperiosteal administration depending on the treatment. In general, deep subcutaneous administration is recommended. A sufficient tissue cover and support are important parameters to achieve a good esthetical treatment result. A minimum of 1 cm of the skin thickness, including subcutaneous fat, is usually required to achieve excellent results and long lasting effect.





Product Features

- 1 graduated sterile pre-filled syringes of 10 ml



8 – 12 months

Treatment

- glutes
- calfs
- body depression



The times indicated should be associated with the patient's lifestyle

Estilad BOOSTER+

A purified, viscous and transparent gel as part of a soft Cross-linked Hyaluronic Acid, the main component which augments facial volume without surgery.

Estilad® BOOSTER+ offers the potential to restore youthful appearance, shaping and tissue fullness with a minimally invasive aesthetic procedure. Estilad® BOOSTER+ is a soft tissue restorative agent and filler, it is primarily used in injection into superficial layer of the subcutaneous tissue on the face to correct minor facial lines and appearance of aged skin. Estilad® BOOSTER+ special formulation is suitable when we want to achieve immediate results with less lasting effect. With special bio-resorbable properties, the filler is ultra-soft.





Product Features

- 1 graduated sterile pre-filled syringes of 1 ml
- 1 x 1 ml and 2 needles 30G



**1 session a week
for 4 times**

Treatment

- anti aging
- damages from photo and
chrono aging

The times indicated should be associated with the patient's lifestyle



Estilad Fine

Pack

1 x 1 ml

Volumizing effect*



PH

7

HA Concentration

21 mg/ml

Needle

2 x 30G

Estilad Touch

Pack

1 x 1 ml

Volumizing effect*



PH

7

HA Concentration

22 mg/ml

Needle

1 x 27G

Estilad Volume

Pack

1 x 1 ml

Volumizing effect*



PH

7

HA Concentration

24 mg/ml

Needle

1 x 27G

*● = light ●● = medium ●●● = strong ●●●● = extra strong ○ = biorevitalising ■ = wide volume

| Estilad Volumax |
|--------------------|
| Pack |
| 1 x 1 ml |
| Volumizing effect* |
| ● ● ● ● |
| PH |
| 7 |
| HA Concentration |
| 26 mg/ml |
| Needle |
| 1 x 26G |

| Estilad Booster+ |
|--------------------|
| Pack |
| 1 x 1 ml |
| Volumizing effect* |
| ○ |
| PH |
| 7 |
| HA Concentration |
| 21 mg/ml |
| Needle |
| 2 x 30G |

| Estilad Contour |
|--------------------|
| Pack |
| 1 x 10 ml |
| Volumizing effect* |
| ■ |
| PH |
| 7 |
| HA Concentration |
| 20 mg/ml |
| Needle |
| - |

*● = light ●● = medium ●●● = strong ●●●● = extra strong ○ = biorevitalising ■ = wide volume



Address

Vancouver

Vancouver, BC V6C 0A3 Canada

Barcelona

Avenida diagonal 539-541, 2izq, Barcelona 08029, Spain

Larnaca

Spetson street 7, 6018, Prodromos, Larnaca - Cyprus

Phone

Spain +34 930 129 369

Cyprus +357 96 694 186

Canada +180 77 882 249

E-mail

office@celiraz-int.net

info@estilad.net

Web site

www.celiraz.com

www.estilad.com

